

NEW

DIGITAL OVERHEAD STIRRER



6000.EU.01

RPM
50-2200 RPM

Max. Capacity
40 Ltr.

Body Material
Aluminum & ABS

Glassco Laboratory Equipments Pvt. Ltd.

CE
certified



info@glasscolabs.com | www.glasscolabs.com
+91 171 2891633, 6451250

50

years of science

DIGITAL OVERHEAD STIRRER *New*


Product Features

- Independent LED display for RPM & Torque real time
- User Friendly
- Accurate display of actual & set value
- Adjustable torque for viscous sample
- Brushless DC motor
- Protects from spillage
- Liquid cannot reach internal components due to front safeguard
- Audible alarm for error conditions
- Effective stirring provides clock wise & anti clock wise stirring modes
- Prominent used in chemical synthesis, pharmaceuticals, physical & chemical analysis, petrochemical, cosmetics, health products, foods, biotechnology & so on

Optional Accessories

Technical Specification :

Body Material	: Aluminum & ABS
Display Window	: RPM & Time
Speed Range	: 50-2200RPM
Speed Control Accuracy	: ± 3 RPM
Speed Display Accuracy	: ± 1 RPM
Max. Viscosity	: 50000 mPas
Max. Torque	: 60NCM
Time Range	: 30~999Min
Max. Stir Capacity (H₂O)	: 40L
Drill Check Dia Range	: 0.5~10mm
Dimension (WxDxH)	: 83 x 220 x 180mm
IP Class	: IP21
Voltage	: 230/115VAC (50/60 Hz)
Power	: 130W



Model No.	V	W	Hz
6000.EU.01	230/115VAC	130 W	50/60Hz