




MICROSLIDES



CAT. NO. 565.303.01

- Manufactured from the best selected glass sheets
- Optically flat glass for microscopic use
- Ground polished edge
- Clean and ready for use
- Size : 76mm x 26mm \pm 1mm
- Refractive Index = 1.5255 \pm .0015

MICROSLIDES FROSTED



CAT. NO. 565.303.02

- Manufactured from the best selected glass sheets
- optically flat glass for microscopic use
- ground polished edge
- clean and ready for use
- Size : 76mm x 26mm \pm 1mm
- Refractive Index = 1.5255 \pm .0015
- Frosted End : 19mm
- Refractive Index = 1.5255 \pm .0015