

NEW

# ROLLER MIXER



RPM  
10-80 RPM

Max. Capacity  
4 kgs

Body Material  
Plastic

Glassco Laboratory Equipments Pvt. Ltd.

CE  
certified



[info@glasscolabs.com](mailto:info@glasscolabs.com) | [www.glasscolabs.com](http://www.glasscolabs.com)  
+91 171 2891633, 6451250

50

years of science


## ROLLER MIXER *New*

### Product Features

- Efficient Rolling and swinging motion for mixing blood samples, Liquid solid suspension & various substance
- Digital set-up for RPM & time with independent LED display
- Maximum stirring speed upto 80RPM
- Seven rollers with small foot print & easy cleaning
- Brushless DC motor
- Compatible for use with wide range of tubes
- Used for prevention of blood coagulation and immune precipitation etc

### Technical Specification :

<b>Mixer Body Material</b>	: Plastic
<b>Operation Mode</b>	: Rolling & Swinging
<b>Amplitude</b>	: 24mm
<b>Roller Quantity</b>	: 7
<b>Speed Range</b>	: 10-80RPM
<b>Speed Unit</b>	: 1RPM
<b>Roller Length</b>	: 280mm
<b>Motor</b>	: brushless DC motor
<b>Dimension (DxWxH)</b>	: 500 x 280 x 120mm
<b>Max. Capacity</b>	: 4kgs
<b>Weight of Body</b>	: 6kg
<b>Voltage</b>	: AC100-240 50-60Hz
<b>Power</b>	: 50W



Model No.	V	W	Hz
4000.EU.01	230 V	50 W	50