

Safety. Performance. Protection that lasts.

3.3
BORO
GLASS

Manufactured from EU-certified
borosilicate 3.3 expansion glass.



External polymer coating
provides protection from
scratches and splinters.



Reduce breakage hazards and
save valuable chemicals.



Safety coating has a high
chemical resistance.



Secure storage of high-potency active
pharmaceutical ingredient (HPAPI)
powders with a safe, reliable grip.



UV amber coating protects against
light up to 450 nm, helping
extend product shelf life.



Lot number (Retrace Code)
for traceability.



LabSafe™ by Glassco

Because safety is not an option.
It's responsibility




SUSTAINABLE BY DESIGN

Coating process is free from
VOCs and isocyanates



Coating helps **prevent
shattering & splashes**



DISCOVERY BEGINS WITH CLARITY



LabSafe™

Safety Coated Bottles

Engineered to protect what matters most.



www.glasscolabs.com





Bottle, LabSafe™ Clear with Screw Cap

Item Code	Capacity (ml)	O.D. x Height (mm)	Screw Thread	Pack of
274.CG.01	100	56x100	GL-45	10
274.CG.02	250	70x138	GL-45	10
274.CG.03	500	86x176	GL-45	10
274.CG.04	1000	101x225	GL-45	10
274.CG.05	2000	136x260	GL-45	10
274.CG.07	5000	181x330	GL-45	1
274.CG.08	10000	227x410	GL-45	1
274.CG.09	20000	288x505	GL-45	1

3.3

Boro Glass



ENQUIRE

*Caps should be loosened while autoclaving



Bottle, LabSafe™ Amber With Screw Cap

Item Code	Capacity (ml)	O.D. x Height (mm)	Screw Thread	Pack of
275.CG.01	100	56x100	GL-45	10
275.CG.02	250	70x138	GL-45	10
275.CG.03	500	86x176	GL-45	10
275.CG.04	1000	101x225	GL-45	10
275.CG.05	2000	136x260	GL-45	10
275.CG.07	5000	183x330	GL-45	1

3.3

Boro Glass



ENQUIRE

*Caps should be loosened while autoclaving



Cap for Bottle GL 45

- Made from Polypropylene
- Autoclavable upto 140° C

Item Code	Colour	Thread	Pack of
275.205.01	Blue	GL 45	10
275.205.01O	Orange	GL 45	10
275.205.01Y	Yellow	GL 45	10
275.205.01G	Green	GL 45	10
275.205.01E	Grey	GL 45	10
275.205.01W	White	GL 45	10
275.205.01B	Black	GL 45	10
275.205.01R	Red	GL 45	10
275.205.00A	Blue	GL 32	10

Pouring Ring for Bottle GL 45

Item Code	Colour	Thread	Pack of
275.205.02	Blue	GL 45	10
275.205.02B	Black	GL 45	10
275.205.02E	Grey	GL 45	10
275.205.02G	Green	GL 45	10
275.205.02O	Orange	GL 45	10
275.205.02R	Red	GL 45	10
275.205.02W	White	GL 45	10
275.205.02Y	Yellow	GL 45	10
275.205.00B	Blue	GL 32	10



ENQUIRE



ENQUIRE

HPLC Safety Cap



Tamper Evident Proof Screw Cap



PTFE Screw Cap



PP Screw Cap




Screw Finish



PTFE Septa

Anti-drip Poring Ring



LabSafe™ Safety Coated



SCAN FOR WEBSITE

GLASSCO INDIA PVT. LTD.

Head Office:
Glassco Laboratory Equipments Pvt. Ltd.
Manglai, P. O. Khuda Kalan, Ambala Cantt. 133 004 (INDIA)

+91 171-2891633 | 70820-15113, 70820-15114

contactus@glasscolabs.com

www.glasscolabs.com