


DEPENDABLE RESULTS
BEGINS WITH
CLARITY



Dispo-LabGlass®

Disposable
Culture Tubes

5.1
Boro
Glass

🌐 www.glasscolabs.com



Manufactured from **Type-I, Class B Borosilicate Glass** for superior clarity, strength and chemical durability.

IDEAL FOR
A WIDE RANGE OF LABORATORY
APPLICATIONS

MICROBIOLOGY & CULTURE

HEMATOLOGY & BLOOD BANKS

CLINICAL & DIAGNOSTIC LABS

ACADEMIC & RESEARCH LABS

INDUSTRIAL QUALITY CONTROL



Chemical
Resistant



High Heat
Tolerance



Convenient
Packaging



Smooth
Finish




Reliable
Performance

Disposable Glassware Culture Tube – Boro Glass 5.1

Made from Type-I borosilicate glass for excellent chemical resistance, heat tolerance, and reliable laboratory use.

- Uniform dimensions & high durability
- Heat-treated above 500°C (932°F)
- Fire-polished tops & rounded bottoms
- Clean-sealed for contamination protection
- 6 x 50 mm size also available as a Durham Tube

US
SIZES



Test Tube Borosilicate 5.1 Expansion Glass

Item Code	O.D. x Ht (mm)	Over Flow Vol (mm)	Pieces/Case
27729.566	6 x 50	0.75	2000
27729.568	10 x 75	3	1000
27729.570	12 x 75	5	1000
27729.572	13 x 100	9	1000
27729.574	15 x 85	10	1000
27729.576	16 x 100	14	1000
27729.578	16 x 125	18	1000
27729.580	16 x 150	22	1000
27729.583	18 x 150	28	500
27729.584	20 x 150	34	500
27729.586	25 x 150	60	500

Ordering Information: Due to machine packaging, quantities per case are approximate.

*Customized sizes also available on request

5.1
Boro
Glass

Chemical
Resistant

White
Printing



ENQUIRE

US
SIZES



With Screw Cap Finish - Boro Glass 5.1

Tubes have a marking spot and reference mark for orientation during use and a threaded neck to accept screw cap.

Salient Features:

- Borosilicate glass 5.1 expansion construction
- Marking spot and orientation reference mark on all tubes

Item Code	O.D. x Ht (mm)	Thread Finish	Over Flow Vol (mm)	Pieces/Case
23283.800	13 x 100	13-415	8	1000
23283.802	16 x 100	15-415	12	1000
23283.804	16 x 125	15-415	16	1000
23283.806	16 x 150	15-415	20	1000
23283.685	20 x 125	18-415	30	500
23283.686	20 x 150	18-415	33	500

*Customized thread finish available in 400, 414 and 425


5.1
Boro
Glass

Chemical
Resistant

White
Printing



ENQUIRE



Phenolic Caps, With Natural Rubber Liner

Item Code	O.D. x Ht (mm)	Over Flow Vol (mm)	Pieces/Case
23283.800C	13 x 100	13-415	1000
23283.804C	16 x 125	15-415	1000
23283.685C	20 x 125	18-415	1000


Ordering Information:

Due to machine packaging, quantities per case are approximate.






ENQUIRE

Disposable Culture Tube (Boro 5.1)



With Rim

Screw Finish



Cork

Phenolic Cap



SCAN FOR WEBSITE

GLASSCO INDIA PVT. LTD.

Head Office:
Glassco Laboratory Equipments Pvt. Ltd.
Manglai, P. O. Khuda Kalan, Ambala Cantt. 133 004 (INDIA)

+91 171-2891633 | 70820-15113, 70820-15114

contactus@glasscolabs.com

www.glasscolabs.com