


# HEATING MANTLE



1000.EU

Temp  
upto 400°C

Mantle Capacity  
50ml - 20L

Heating Mantle Material  
Glass Fibre Yarn

Glassco Laboratory Equipments Pvt. Ltd.

CE  
certified



[info@glasscolabs.com](mailto:info@glasscolabs.com) | [www.glasscolabs.com](http://www.glasscolabs.com)  
+91 171 2891633, 6451250

45

years of science



# HEATING MANTLE

- Aluminium Moulded Construction
- ON-OFF knob for heat power control.
- Mantle knitted with glass fibre yarn.
- Provides uniform heating to flasks.
- Heats upto 400°C.
- For Round Bottomed Flasks as per DIN - 12347
- Inbuilt protection from liquid spill (upto 5000ml)
- Easy maintenance
- Long life



Liquid Spill Protection

## Technical Specifications

Heating Mantle Body	: Made of aluminium
Heating Mantle Material	: Glass Fibre Yarn
Mantle Capacity	: To hold standard Round Bottom Flasks upto 20 liters
Temperature Range	: Approx. 400°C
Heating Control	: Electronic power regulator
Speed Control	: Electronic regulator/ Energy regulator
Permissible Ambient Temperature	: 5- 40 °C
Permissible relative moisture	: 80%
IP Class	: 10

Model No. (230 VAC 50Hz)	Model No. (115 VAC 60Hz)	Mantle Capacity (ml)	Wattage (W)
1000.EU.01	1000.US.01	50ml	150
1000.EU.02	1000.US.02	100ml	150
1000.EU.03	1000.US.03	250ml	150
1000.EU.04	1000.US.04	500ml	200
1000.EU.05	1000.US.05	1000ml	300
1000.EU.06	1000.US.06	2000ml	450
1000.EU.07	1000.US.07	3000ml	450
1000.EU.08	1000.US.08	4000ml	450
1000.EU.09	1000.US.09	5000ml	600
1000.EU.10	1000.US.10	6000ml	600
1000.EU.11	1000.US.11	10000ml	1200
1000.EU.12	1000.US.12	20000ml	1800